

# KLM Orthotic Labs System\_Rx Sizing Chart - Please use the System\_Rx Order Form to place orders

(Relative) Shoe Size	M1	M2	M3	M4	M5-6	M7-8	M9	M10-11	M12-13	M14	M15-16
Length	W3	W4	W5	W6	W7-8	W9-10	W11	W13-14	W14-15	W16	W17-18
	118	123	128	138	148	158	168	178	188	198	208
	<b>Heel 40</b> <b>X019 (1)</b> 118-40-3-14	<b>Heel 42</b> <b>X020 (2)</b> 123-42-3-14	<b>Heel 38</b> <b>X126</b> 128-38-3-17	<b>Heel 47</b> <b>X022 (4)</b> 138-47-3-14	<b>Heel 47</b> <b>X014</b> 148-47-3-14	<b>Heel 52</b> <b>X005</b> 158-52-3-14	<b>Heel 52</b> <b>X008</b> 168-52-3-14	<b>Heel 52</b> <b>X018</b> 178-52-3-14	<b>Heel 57</b> <b>X013</b> 188-57-3-14	<b>Heel 62</b> <b>X015</b> 198-62-3-14	<b>Heel 67</b> <b>X131</b> 208-67-3-17
		<b>X129</b> 123-42-3-17	<b>Heel 42</b> <b>X026</b> 128-42-3-17	<b>X123</b> 138-47-3-17	<b>X118 (5)</b> 148-47-3-17	<b>X101 (7)</b> 158-52-3-17	<b>X108</b> 168-52-3-17	<b>X117</b> 178-52-3-17	<b>X113</b> 188-57-3-17	<b>X115</b> 198-62-3-17	<b>X222</b> 208-67-3-20
		<b>Heel 45</b> <b>X025</b> 123-45-3-14	<b>X127</b> 128-42-3-17	<b>Heel 52</b> <b>X024</b> 138-52-3-14	<b>Heel 52</b> <b>X009</b> 148-52-3-14	<b>X201</b> 158-52-3-20	<b>X208</b> 168-52-3-20	<b>X215</b> 178-52-3-20	<b>X214</b> 188-57-3-20	<b>X212</b> 198-62-3-20	<b>X319</b> 208-67-3-23
		<b>X130</b> 123-45-3-17	<b>Heel 45</b> <b>X021 (3)</b> 128-45-3-14	<b>X120</b> 138-52-3-17	<b>X111 (6)</b> 148-52-3-17	<b>X313</b> 158-52-3-23	<b>X314</b> 168-52-3-23	<b>X315</b> 178-52-3-23	<b>X310</b> 188-57-3-23	<b>X309</b> 198-62-3-23	
			<b>X122</b> 128-45-3-17	<b>X216</b> 148-52-3-20	<b>X002</b> 158-57-3-14	<b>Heel 57</b> <b>X003</b> 168-57-3-14	<b>Heel 57</b> <b>X007</b> 178-57-3-14	<b>Heel 57</b> <b>X006</b> 188-62-3-14	<b>Heel 62</b> <b>X007</b> 198-67-3-14	<b>Heel 67</b> <b>X016</b> 208-72-3-17	<b>Heel 72</b> <b>X125</b> 208-72-3-23
			<b>Heel 47</b> <b>X023</b> 128-47-3-14	<b>X017</b> 148-57-3-14	<b>X105 (8)</b> 158-57-3-17	<b>X102</b> 168-57-3-17	<b>X106</b> 178-57-3-17	<b>X107</b> 188-62-3-17	<b>X107</b> 198-67-3-17	<b>X114</b> 208-72-3-20	<b>X219</b> 208-72-3-23
			<b>X119</b> 128-47-3-17	<b>X116</b> 148-57-3-17	<b>X205</b> 158-57-3-20	<b>X202</b> 168-57-3-20	<b>X206</b> 178-57-3-20	<b>X207</b> 188-62-3-20	<b>X207</b> 198-67-3-20	<b>X211</b> 208-72-3-23	<b>X318 (15)</b> 208-72-3-23
				<b>X217</b> 148-57-3-20	<b>X311</b> 158-57-3-23	<b>X301</b> 168-57-3-23	<b>X306</b> 178-57-3-23	<b>X304</b> 188-62-3-23	<b>X304</b> 198-67-3-23	<b>X308</b> 208-72-3-23	
					<b>Heel 62</b> <b>X011</b> 158-62-3-14	<b>Heel 62</b> <b>X001</b> 168-62-3-14	<b>Heel 62</b> <b>X004</b> 178-62-3-14	<b>Heel 62</b> <b>X010</b> 188-67-3-14	<b>Heel 67</b> <b>X010</b> 198-72-3-17	<b>Heel 72</b> <b>X121</b> 208-77-3-17	<b>Heel 77</b> <b>X128</b> 208-77-3-20
					<b>X112</b> 158-62-3-17	<b>X104 (9)</b> 168-62-3-17	<b>X103 (10)</b> 178-62-3-17	<b>X109</b> 188-67-3-17	<b>X109</b> 198-72-3-20	<b>X220</b> 208-77-3-20	<b>X221</b> 208-77-3-23
					<b>X213</b> 158-62-3-20	<b>X204</b> 168-62-3-20	<b>X203</b> 178-62-3-20	<b>X209 (12)</b> 188-67-3-20	<b>X209 (12)</b> 198-72-3-23	<b>X316 (14)</b> 208-77-3-23	<b>X320 (16)</b> 208-77-3-23
					<b>X312</b> 158-62-3-23	<b>X303</b> 168-62-3-23	<b>X302</b> 178-62-3-23	<b>X305</b> 188-67-3-23	<b>X305</b> 198-72-3-23		
							<b>Heel 67</b> <b>X012</b> 178-67-3-14	<b>Heel 72</b> <b>X124</b> 188-72-3-17	<b>Heel 72</b> <b>X218</b> 188-72-3-20		
							<b>X110</b> 178-67-3-17	<b>X218</b> 188-72-3-20	<b>X317 (13)</b> 188-72-3-23		
							<b>X210 (11)</b> 178-67-3-20				
							<b>X307</b> 178-67-3-23				



For more information contact us at:  
 phone: 800-556-3668  
 fax: 800-556-3338  
 email: cservice@klmlabs.com

- ### Sizing Procedure
- 1) Determine length by shoe size
  - 2) Determine heel width
  - 3) Determine medial arch height

\* All measurements are in millimeters

- ### Choose Shell Stiffness/Color
- Standard Rx Shell (Semi-Rigid) Blue, Black, White
  - Rx Rigid-Flex (In Between Semi-Rigid and Semi\_Flexible) Brown
  - Rx Semi-Flex (Semi-Flexible) Flesh Tone
  - Rx Flex (Flexible) Clear

